
Rescue Manual

Zeramex XT Rescue Set

The Zeramex XT Rescue Set is a set of instruments designed to remove a broken or worn abutment screw.

Table of contents

Rescue Set – Product Overview 2

Case 1: Removal of a worn Vicarbo screw through the abutment..... 3

Case 2: Removal of a fractured Vicarbo screw..... 5

General information 7



 Precautions: 7


 Sterilization..... 7

 Removing the screw from the Rescue Instrument 7

 Rescue instrument check prior to use 7

Rescue Set – Product Overview


	Drilling pilot hole				Screw removal		Thread cleaner
							
Included in Rescue Set	yes	yes	yes	Please note: Not included in the Rescue Set	yes	Please note: Not included in the Rescue Set	Please note: Not included in the Rescue Set
Article No.	R98120	R98100	R98110	*Pilot Drill	R98130	P48935	Interdental brush
Name	Rescue Shell XT Small	Rescue Shell XT Regular	Rescue Shell XT Wide	*Pilot Drill	Rescue Instrument 28mm	Surgical Ratchet	Interdental brush (e.g. from Peridenta Care)
Material	PEEK	PEEK	PEEK	n/a	Stainless steel	Stainless steel	Nylon brush

 The electronic Instructions for Use for our products are available at ifu.zeramex.com


*Pilot Drill - Not included in the Rescue-Set. Specification:

Item: Standard Carbide Bur/Drill
Tip Type: Round burr, \varnothing 1,0 mm (O10)
Shank Type: Right Angle Long (RAL) or (RAXL)
Shank: RAL (or longer)

e.g. Komet item No. [H1SEM.205.010](#) / Meisinger USA - [HM1S-023-RAL](#)

Case 1: Removal of a worn Vicarbo screw through the abutment

Depth marking of the drill

Provide a pilot drill according to the specification on page 2.

[1] Guide the Pilot drill into the Abutment until it reaches the end position and mark [2] the shaft of the drill with a waterproof pen.

Remove the drill bit and mark a line 3-3.5mm above the first mark [3], to visualize the final position of the drilling depth.

Hint: The drilling depth corresponds to the Rescue Instrument thread length.

Drill (Pilot hole 3.5mm)

Insert the Pilot drill into the handpiece.


[4] Guide the Pilot drill through the abutment. Drill a 3-3.5mm deep pilot hole exactly in the center of the screw at a speed of 1,000 rpm intermittently and with as little pressure as possible.

The drilling process is complete when the marking [5] is flush with the top of the abutment chimney.

CAVE: Always use a new pilot drill for each pilot hole.

Note: Use a dental dam if suctioning the particles is not sufficient.


🔄 opt. 1'000 rpm

Pilot drill (e.g. Komet item No. [H1SEM.205.010](#)
Meisinger USA - [HM1S-023-RAL](#))


Tip: The drilling process is simplified by applying slight pressure and pulling the drill bit out again. Ensure centric drilling.

Case 1: Removal of a worn Vicarbo screw through the abutment


Screw removal

[6] Guide the Rescue Instrument into the pilot hole of the screw and turn the adapter by hand or ratchet [7] to the left (counterclockwise) while maintaining pressure until the screw is removed [8].

Cleaning

Clean the interior of the implant with the cleaning brush [9] and handpiece. Make sure that there are no black residues in the implant. A rotating interdental brush (e.g. from Peridenta Care¹) with the contra-angle handpiece is recommended for cleaning.


If the screw cannot be removed

1. Check the Rescue Instrument thread for signs of wear as shown in the illustration in the last chapter.
2. If the Rescue Instrument fails to engage with the screw-/fragment repeat the drilling step while making sure that the pilot hole is centric and sufficiently deep. The pilot hole can also be drilled without the Rescue Shell if it is considered to be an obstacle and you have a clear view of the fragment.
3. If the fragment still cannot be removed by the Rescue Instrument, the pilot drill can be used to drill out the entire fragment as a last resort. Make sure that adequate cooling is provided.

¹ Not included in set.

Case 2: Removal of a fractured Vicarbo screw

Depth marking of the drill


Provide a pilot drill according to the specification on page 2.

[1] Place the Rescue Shell Small, Regular or Wide on the corresponding implant and make sure it is seated securely. A light-cured composite (e.g. primo pattern LC) can be used to fix the Rescue Shell in connection with the adjacent teeth.

[2] Guide the Pilot drill into the Rescue shell until it reaches the end position and [3] mark the shaft of the drill with a waterproof pen.

[4] Remove the drill bit and mark a line 3-3.5mm above the first mark, to visualize the final position of the drilling depth. Hint: The drilling depth corresponds to the Rescue Instrument thread length.

Note: If the guide sleeve is perceived as an obstacle and the circumstances and circumstances allow a good view of the screw fragment, the centric pilot hole can also be drilled freehand.

Pilot drill (e.g. Komet item No. [H1SEM.205.010](#) Meisinger USA - [HM1S-023-RAL](#))

Drill (Pilot hole 3.5mm)

Insert the Pilot drill into the handpiece.

[5] Guide the Pilot drill through the Rescue Shell.


Drill a 3-3.5mm deep pilot hole exactly in the center of the screw at a speed of 1,000 rpm intermittently and with as little pressure as possible.

The drilling process is complete when the marking [6] is flush with the top plane of the Rescue shell.

Note: Use a dental dam if suctioning the particles is not sufficient.


CAVE: Always use a new pilot drill for each pilot hole.

opt. 1'000 rpm


Tip: The drilling process is simplified by applying slight pressure and pulling the drill bit out again. Ensure centric drilling.

Case 2: Removal of a fractured Vicarbo screw


Screw removal

[7] Guide the Rescue Instrument into the pilot hole of the screw and turn the adapter by hand or ratchet [8] to the left (counterclockwise) while maintaining pressure until the screw is removed [9].

Cleaning

[10] Clean the interior of the implant with the cleaning brush and handpiece. Make sure that there are no black residues in the implant. A rotating interdental brush (e.g. from Peridenta Care²) with the contra-angle handpiece is recommended for cleaning.

If the screw cannot be removed

1. Check the Rescue Instrument thread for signs of wear as shown in the illustration in the last chapter.
2. If the Rescue Instrument fails to engage with the screw-/fragment repeat the drilling step while making sure that the pilot hole is centric and sufficiently deep. The pilot hole can also be drilled without the Rescue Shell if it is considered to be an obstacle and you have a clear view of the fragment.
3. If the fragment still cannot be removed by the Rescue Instrument, the pilot drill can be used to drill out the entire fragment as a last resort. Make sure that adequate cooling is provided.

² Not included in set.

General information

Precautions:

- Our products must be protected against aspiration when handled intraorally.
- Non-sterile packaged products must not be sterilized in their original packaging.

Sterilization

Devices delivered non-sterile must be cleaned and disinfected prior to use. Prior to processing, remove the packaging system from the product.


Devices delivered non-sterile must be sterilized prior to use. For further details on sterilization procedures, please consult the IFU online at ifu.zeramex.com

Removing the screw from the Rescue Instrument

To remove the screw from the Rescue Instrument, hold the screw in the threaded area with pliers and turn it clockwise with the ratchet.

Rescue instrument check prior to use

Check the thread for wear, especially on the lower half of the thread tip. Rescue instrument thread flanks with signs of wear, such as those on the image below (marked with a red cross), must not be used. Otherwise, the function is impaired, and damage may occur.

Symbols

The following table describes the symbols that may be printed on the packaging label. Please refer to the packaging label for the applicable symbols related to the product.

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