

# Importing CAD libraries Zeramex XT

## How to import the library into your **exocad** system?

1. Close your exocad software
2. Unzip the folder “zeramex\_xt\_exocad\_v1.zip”
3. Move the unzipped folder to the following installation directory of the exocad software: \DentalCADApp\library\implant
4. Launch your exocad software
5. You have successfully installed the Zeramex XT Zerabase library.

**Information:** Zeramex XT exocad Model Creator will be available soon.

## How to import the library into your **3Shape** system?

1. Open the 3Shape Dental System Control Panel
2. Tools ► Import/Export ► Import material
3. Select “zeramex\_xt\_3shape\_v5.dme”. Save and close the “3Shape Dental System Control Panel”.
4. You have successfully installed the Zeramex XT Zerabase library.

Please **contact the support service of your CAD/CAM provider with any questions** on software data import.

Order the parts necessary for treatment directly with Zeramex.

### Order-Hotline

T 00800 – 93 55 66 37

F 00800 – 93 55 63 77

[order@zeramex.com](mailto:order@zeramex.com)