






ZERAMEX<sup>®</sup> (P6)



Sortiment

## ZERAMEX® P6 Implants


Ø3.3mm	Art.Nr.	Name	Dimension	Material
	P15508	ZERAMEX®P6 Implant Ø3.3x8mm SN (incl. Healing Cap)	Länge: 8mm	ZrO <sub>2</sub> -ATZ-HIP
	P15510	ZERAMEX®P6 Implant Ø3.3x10mm SN (incl. Healing Cap)	Länge: 10mm	
	P15512	ZERAMEX®P6 Implant Ø3.3x12mm SN (incl. Healing Cap)	Länge: 12mm	
Ø4.1mm	Art.Nr.	Name	Dimension	Material
	P16508	ZERAMEX®P6 Implant Ø4.1x8mm RN (incl. Healing Cap)	Länge: 8mm	ZrO <sub>2</sub> -ATZ-HIP
	P16510	ZERAMEX®P6 Implant Ø4.1x10mm RN (incl. Healing Cap)	Länge: 10mm	
	P16512	ZERAMEX®P6 Implant Ø4.1x12mm RN (incl. Healing Cap)	Länge: 12mm	
Ø4.8mm	Art.Nr.	Name	Dimension	Material
	P17508	ZERAMEX®P6 Implant Ø4.8x8mm RN (incl. Healing Cap)	Länge: 8mm	ZrO <sub>2</sub> -ATZ-HIP
	P17510	ZERAMEX®P6 Implant Ø4.8x10mm RN (incl. Healing Cap)	Länge: 10mm	
	P17512	ZERAMEX®P6 Implant Ø4.8x12mm RN (incl. Healing Cap)	Länge: 12mm	





## ZERAMEX® P6 Abutments


SN Plattform	Art.Nr.	Name	Dimension	Material
	P15501	ZERAMEX®P6 Abutment SN Straight, 1mm (incl. Screw)	Ø: 5.0mm, AH: 6.5mm, KH: 1mm	ZrO <sub>2</sub> -ATZ-HIP VICARBO®
	P15502	ZERAMEX®P6 Abutment SN Straight, 2mm (incl. Screw)	Ø: 5.0mm, AH: 7.5mm, KH: 2mm	
	P15515	ZERAMEX®P6 Abutment SN Angular 15°, 1mm (incl. Screw)	Ø: 5.0mm, AH: 5.5mm	ZrO <sub>2</sub> -ATZ-HIP VICARBO®

### CAD/CAM

	P15530	ZERAMEX®P6 ZERABASE SN, for crown (incl. Screw)	Ø: 5.0mm, AH: 7mm, KH: 1mm	ZrO <sub>2</sub> -ATZ-HIP VICARBO®
	P15531	ZERAMEX®P6 ZERABASE SN, for bridge (incl. Screw)	Ø: 5.0mm, AH: 7mm, KH: 1mm	


RN Plattform	Art.Nr.	Name	Dimension	Material
	P16501	ZERAMEX®P6 Abutment RN Straight, 1mm (incl. Screw)	Ø: 5.8mm, AH: 7mm, KH: 1mm	ZrO <sub>2</sub> -ATZ-HIP VICARBO®
	P16502	ZERAMEX®P6 Abutment RN Straight, 2mm (incl. Screw)	Ø: 5.8mm, AH: 8mm, KH: 2mm	
	P16515	ZERAMEX®P6 Abutment RN Angular 15°, 1mm (incl. Screw)	Ø: 5.8mm, AH: 7mm, KH: 1mm	ZrO <sub>2</sub> -ATZ-HIP VICARBO®

### CAD/CAM

	P16530	ZERAMEX®P6 ZERABASE RN, for crown (incl. Screw)	Ø: 5.6mm, AH: 7mm, KH: 1mm	ZrO <sub>2</sub> -ATZ-HIP VICARBO®
	P16531	ZERAMEX®P6 ZERABASE RN, for bridge (incl. Screw)	Ø: 5.6mm, AH: 7mm, KH: 1mm	





ZERABASE Abutments sind in die Systeme von 3Shape und Exocad integriert

VICARBO® Schrauben	Art.Nr.	Name	Dimension	Material
	P15450	ZERAMEX®P6 VICARBO SN Straight	Länge: 7.3mm	VICARBO®
	P15451	ZERAMEX®P6 VICARBO SN Angular	Länge: 6.25mm	
	P16450	ZERAMEX®P6 VICARBO RN Straight	Länge: 7.5mm	
	P16451	ZERAMEX®P6 VICARBO RN Angular	Länge: 6.7mm	

## ZERAMEX® ZERALOCK™ Abutments

### Gerade & abgewinkelte Abutments

SN Plattform	Art.Nr.	Name	Dimension	Material
	P15305	ZERAMEX® ZERALOCK™ Abutment SN straight	AH: 6mm	ZrO <sub>2</sub> -ATZ-HIP
	T15315	ZERAMEX® ZERALOCK™ Abutment SN angular 15°	AH: 5.4mm	

RN Plattform	Art.Nr.	Name	Dimension	Material
	P16305	ZERAMEX® ZERALOCK™ Abutment RN straight	AH: 6mm	ZrO <sub>2</sub> -ATZ-HIP
	T16315	ZERAMEX® ZERALOCK™ Abutment RN angular 15°	AH: 5.5mm	





### LOCATOR®



SN Plattform	Art.Nr.	Name	Dimension	Material
	P15202	ZERAMEX® ZERALOCK™ Locator® 2mm SN	Höhe: 2mm	ZrO <sub>2</sub> -ATZ-HIP
	P15203	ZERAMEX® ZERALOCK™ Locator® 3mm SN	Höhe: 3mm	
	P15205	ZERAMEX® ZERALOCK™ Locator® 5mm SN	Höhe: 5mm	

RN Plattform	Art.Nr.	Name	Dimension	Material
	P16202	ZERAMEX® ZERALOCK™ Locator® 2mm RN	Höhe: 2mm	ZrO <sub>2</sub> -ATZ-HIP
	P16203	ZERAMEX® ZERALOCK™ Locator® 3mm RN	Höhe: 3mm	
	P16205	ZERAMEX® ZERALOCK™ Locator® 5mm RN	Höhe: 5mm	






## ZERAMEX® Prosthetics - SN Plattform

Weichgewebemanagement	Art.Nr.	Name	Dimension	Material
	P35500	ZERAMEX®P6 Healing Cap SN	Höhe: 0.5mm	PEEK
	P35503	ZERAMEX®P6 Gingivaformer SN, 3mm	Höhe: 3mm	
	P35504	ZERAMEX®P6 Gingivaformer SN, 4mm	Höhe: 4mm	
	P35530	ZERAMEX®P6 Provisorium SN (incl. Screw)	AH: 7mm KH: 1mm Ø: 4.5mm	PEEK, PEEK-CW30


Abdrucknahme	Art.Nr.	Name	Dimension	Material
	P35510	ZERAMEX®P6 Transfer Open Tray SN	Höhe Hülse: 10mm H. inkl. Schraube: 19mm	PEEK-CW30, Aluminium
	P35512	ZERAMEX®P6 Transfer Closed Tray SN	Höhe Hülse: 7mm H. inkl. Schraube: 14mm	
	P35513	ZERAMEX®P6 Transfer Closed Tray, long SN	Höhe Hülse: 11mm H. inkl. Schraube: 19mm	
	P35514	ZERAMEX®P6 Scanbody SN	Höhe: 10mm	PEEK, PEEK-CW30




# SORTIMENT





## ZERAMEX® Prosthetics - SN Plattform

Hilfsteile	Art.Nr.	Name	Dimension	Material
	P35520	ZERAMEX®P6 Implant-Replica SN	Länge: 10mm	Aluminium
	P35550	ZERAMEX®P6 Provisorium Screw SN	Länge: 7.3mm	PEEK-CW30
	P35450	ZERAMEX®P6 Lab Screw Straight SN	Länge: 7.4mm	PEEK
	P35451	ZERAMEX®P6 Lab Screw Angular SN	Länge: 6.3mm	PEEK


## ZERAMEX® Prosthetics - RN Plattform

Weichgewebemanagement	Art.Nr.	Name	Dimension	Material
	P36500	ZERAMEX®P6 Healing Cap RN	Höhe: 0.6mm	PEEK
	P36502	ZERAMEX®P6 Closure Screw RN	Höhe: 0.5mm	
	P36503	ZERAMEX®P6 Gingivaformer RN, 3mm	Höhe: 3mm	
	P36504	ZERAMEX®P6 Gingivaformer RN, 4mm	Höhe: 4mm	
	P36530	ZERAMEX®P6 Provisorium RN (incl. Screw)	AH: 7mm KH: 1mm Ø: 5.5mm	PEEK, PEEK-CW30









Abdrucknahme	Art.Nr.	Name	Dimension	Material
	P36510	ZERAMEX®P6 Transfer Open Tray RN	Höhe Hülse: 11mm H. inkl. Schraube: 20mm	PEEK-CW30, Aluminium
	P36512	ZERAMEX®P6 Transfer Closed Tray RN	Höhe Hülse: 7mm H. inkl. Schraube: 14mm	PEEK-CW30, Aluminium
	P36513	ZERAMEX®P6 Transfer Closed Tray, long RN	Höhe Hülse: 11mm H. inkl. Schraube: 18mm	PEEK-CW30, Aluminium
	P36514	ZERAMEX®P6 Scanbody RN	Höhe: 10mm	PEEK, PEEK-CW30

Hilfsteile	Art.Nr.	Name	Dimension	Material
	P36520	ZERAMEX®P6 Implant-Replica RN	Länge: 10mm	Aluminium
	P36550	ZERAMEX®P6 Provisorium Screw RN	Länge: 6.1mm	PEEK-CW30
	P36450	ZERAMEX®P6 Lab Screw Straight RN	Länge: 7.5mm	PEEK
	P36451	ZERAMEX®P6 Lab Screw Angular RN	Länge: 6.7mm	PEEK
	P18550	ZERAMEX®P6 Abutment Planning Kit		PEEK














## ZERAMEX® Tools



















	Art.Nr.	Name	Material
	P35601	ZERAMEX®P Rosedrill ø2,2mm	rostfreier Stahl
	P35602	ZERAMEX®P Pilot Drill ø2,2mm	
	P35633	ZERAMEX®P Drill ø2,8mm	
	P36633	ZERAMEX®P Drill ø3,5mm	
	P37633	ZERAMEX®P Drill ø4,2mm	
	P35622	ZERAMEX®P Profile Drill SN ø2.8mm	
	P36622	ZERAMEX®P Profile Drill ø3,5mm	
	P35651	ZERAMEX®P Depth Gauge ø2,2mm	
	P36651	ZERAMEX®P Depth Gauge ø3,5mm	
	P35620	ZERAMEX®P Screw Tap ø3,3mm	
	P36620	ZERAMEX®P Screw Tap ø4,1mm	
	P37620	ZERAMEX®P Screw Tap ø4,8mm	
	P38619	ZERAMEX® Prothetikschlüssel, 19mm	
	P38623	ZERAMEX® Prothetikschlüssel, 23mm	
	T38620	ZERAMEX® Pickup 20mm	
	T38622	ZERAMEX® Rescue Pickup	
	T38625	ZERAMEX® Pickup 25mm	
	T38626	ZERAMEX® Ausrichthilfe	
	C7650	Drill extension	

## ZERAMEX® Tools

	Art.Nr.	Name	Material
	P35453	ZERAFIX™ Platzhalter SN 3.0mm	PEEK
	P36453	ZERAFIX™ Platzhalter RN 3.3mm	
	P38645S	ZERALOCK™ Abutmenthalter-Set SN & RN	PEEK
	P38644S	ZERAFIX™ Abutmenthalter-Set SN & RN	
	P48850	ZERAMEX®P Surgery Tray	rostfreier Stahl
	P48854	ZERAMEX®P Surgery Tray, fully equipped with Ratchet	
	P48870	ZERAMEX®P Surgery Tray "mini"	
	P48932	ZERAMEX®P Ratchet Adapter Unit Short	
	P48935	ZERAMEX®P Surgical Ratchet	
	P48860	ZERAMEX®P Prosthetics Tray	
	P48865	ZERAMEX®P Prosthetics Kit, incl. Ratchet	

## ZERAMEX® Locator® Zubehör

	Art.Nr.	Name	Material
	T38219	LOCATOR® Processing Denture Cap (up to 10° divergence per implant, clear, pink, blue), 2pcs.	Ti/Nylon/LDPE
	T38220	LOCATOR® Processing Denture Cap (up to 20° divergence per implant, green, orange, red), 2pcs.	Ti/Nylon/LDPE
	T38222	LOCATOR® Blockout Spacer, 20pcs.	Silikon-Kautschuk
	T38224	LOCATOR® Impression Coping, 4pcs.	Al-Gehäuse mit LDPE-Einsatz
	T38226	LOCATOR® Replica Ø4mm, 4pcs.	Aluminium
	T38227	LOCATOR® Replica Ø5mm, 4pcs.	Aluminium
	T38228	LOCATOR® Core Tool	rostfreier Stahl
	T38230	LOCATOR® Replacement clear, 4pcs. Retentionskraft: 2.27kg	Nylon
	T38231	LOCATOR® Replacement pink, 4pcs. Retentionskraft: 1.36kg	
	T38232	LOCATOR® Replacement blue, 4pcs. Retentionskraft: 0.68kg	
	T38233	LOCATOR® Replacement green, 4pcs. Retentionskraft: 1.81kg	
	T38234	LOCATOR® Replacement orange, 4pcs. Retentionskraft: 0.91kg	
	T38235	LOCATOR® Replacement red, 4pcs. Retentionskraft: 0.45kg	

CE	0050		Vorsicht
	Bestellnummer		Vor dem Verfallsdatum verwenden
	Chargencode		Gebrauchsanweisung beachten
	Nicht wiederverwenden		Herstellungsdatum
	Dampfsterilisiert		Hersteller
	Unsteril		Bei beschädigter Verpackung nicht verwenden
	Vor direkter Sonneneinstrahlung schützen		Die Produkte des Portfolios sind von der FDA zugelassen
	Trocken lagern		«Gemäss US- Bundesgesetz: Verkauf nur direkt an ausgebildete Ärzte, Zahnärzte oder lizenzierte Fachleute bzw. in deren Auftrags»
	Implantat small ø3,3mm SN (Plattform ø4,0mm)		Implantat small ø4,1mm RN (Plattform ø4,8mm)
	Implantat small ø4,8mm RN (Plattform ø4,8mm)		