



ZERAMEX[®] (P6)


Product range

PRODUCT RANGE

ZERAMEX®P6 Implants



Ø 3.3 mm	Art. No.	Name	Dimension	Material
	P15508	ZERAMEX®P6 Implant Ø 3.3 x 8 mm SN (incl. healing cap)	Length: 8mm	ZrO ₂ -ATZ-HIP
	P15510	ZERAMEX®P6 Implant Ø 3.3 x 10 mm SN (incl. healing cap)	Length: 10mm	
	P15512	ZERAMEX®P6 Implant Ø 3.3 x 12 mm SN (incl. healing cap)	Length: 12mm	

Ø 4.1 mm	Art. No.	Name	Dimension	Material
	P16508	ZERAMEX®P6 Implant Ø 4.1 x 8 mm RN (incl. healing cap)	Length: 8mm	ZrO ₂ -ATZ-HIP
	P16510	ZERAMEX®P6 Implant Ø 4.1 x 10 mm RN (incl. healing cap)	Length: 10mm	
	P16512	ZERAMEX®P6 Implant Ø 4.1 x 12 mm RN (incl. healing cap)	Length: 12mm	


Ø 4.8 mm	Art. No.	Name	Dimension	Material
	P17508	ZERAMEX®P6 Implant Ø 4.8 x 8 mm RN (incl. healing cap)	Length: 8mm	ZrO ₂ -ATZ-HIP
	P17510	ZERAMEX®P6 Implant Ø 4.8 x 10 mm RN (incl. healing cap)	Length: 10mm	
	P17512	ZERAMEX®P6 Implant Ø 4.8 x 12 mm RN (incl. healing cap)	Length: 12mm	





ZERAMEX® P6 Abutments

SN Platform	Art. No.	Name	Dimension	Material
	P15501	ZERAMEX®P6 Abutment SN Straight, 1 mm (incl. screw)	AH: 6.5mm, CH: 1mm	ZrO ₂ -ATZ-HIP VICARBO®
	P15502	ZERAMEX®P6 Abutment SN Straight, 2 mm (incl. screw)	AH: 6.5mm, CH: 2mm	
	P15515	ZERAMEX®P6 Abutment SN Angular 15°, 1 mm (incl. screw)	AH: 5.5mm	ZrO ₂ -ATZ-HIP VICARBO®

CAD/CAM


	P15530	ZERAMEX®P6 ZERABASE SN, for crown (incl. screw)	AH: 7mm, CH: 1mm	ZrO ₂ -ATZ-HIP VICARBO®
	P15531	ZERAMEX®P6 ZERABASE SN, for bridge (incl. screw)	AH: 7mm, CH: 1mm	

RN Platform	Art. No.	Name	Dimension	Material
	P16501	ZERAMEX®P6 Abutment RN Straight, 1 mm (incl. screw)	AH: 7mm, CH: 1mm	ZrO ₂ -ATZ-HIP VICARBO®
	P16502	ZERAMEX®P6 Abutment RN Straight, 2 mm (incl. screw)	AH: 8mm, CH: 2mm	
	P16515	ZERAMEX®P6 Abutment RN Angular 15°, 1 mm (incl. screw)	AH: 7mm, CH: 1mm	ZrO ₂ -ATZ-HIP VICARBO®




The ZERABASE Abutments are integrated in the 3Shape and Exocad systems.


AH: Abutment height CH: Collar height

VICARBO® screws	Art. No.	Name	Dimension	Material
	P15450	ZERAMEX®P6 VICARBO SN Straight	Length: 7.3mm	VICARBO®
	P15451	ZERAMEX®P6 VICARBO SN Angular	Length: 6.25mm	
	P16450	ZERAMEX®P6 VICARBO RN Straight	Length: 7.5mm	
	P16451	ZERAMEX®P6 VICARBO RN Angular	Length: 6.7mm	

ZERAMEX® ZERALOCK™ Abutments



Straight & angular abutments

SN Platform	Art. No.	Name	Dimension	Material
	P15305	ZERAMEX® ZERALOCK™ Abutment SN straight	AH: 6mm	ZrO ₂ -ATZ-HIP
	T15315	ZERAMEX® ZERALOCK™ Abutment SN angular 15°	AH: 5.4mm	

RN Platform	Art. No.	Name	Dimension	Material
	P16305	ZERAMEX® ZERALOCK™ Abutment RN straight	AH: 6mm	ZrO ₂ -ATZ-HIP
	T16315	ZERAMEX® ZERALOCK™ Abutment RN angular 15°	AH: 5.5mm	





LOCATOR®

SN Platform	Art. No.	Name	Dimension	Material
	P15202	ZERAMEX® ZERALOCK™ Locator® 2 mm SN	Height: 2mm	ZrO ₂ -ATZ-HIP
	P15203	ZERAMEX® ZERALOCK™ Locator® 3 mm SN	Height: 3mm	
	P15205	ZERAMEX® ZERALOCK™ Locator® 5 mm SN	Height: 5mm	
RN Platform	Art. No.	Name	Dimension	Material
	P16202	ZERAMEX® ZERALOCK™ Locator® 2 mm RN	Height: 2mm	ZrO ₂ -ATZ-HIP
	P16203	ZERAMEX® ZERALOCK™ Locator® 3 mm RN	Height: 3mm	
	P16205	ZERAMEX® ZERALOCK™ Locator® 5 mm RN	Height: 5mm	






ZERAMEX® Prosthetics - SN Platform

Soft Tissue Management	Art. No.	Name	Dimension	Material
	P35500	ZERAMEX®P6 Healing Cap SN	Height: 0.5mm	PEEK
	P35503	ZERAMEX®P6 Gingiva Former SN, 3 mm	Height: 3mm	
	P35504	ZERAMEX®P6 Gingiva Former SN, 4 mm	Height: 4mm	
	P35530	ZERAMEX®P6 Provisional SN (incl. screw)	AH: 7mm CH: 1mm	PEEK, PEEK-CW30

Impression Taking	Art. No.	Name	Dimension	Material
	P35510	ZERAMEX®P6 Transfer Open Tray SN	Sleeve height: 10mm Height incl. screw: 19mm	PEEK-CW30, Aluminum
	P35512	ZERAMEX®P6 Transfer Closed Tray SN	Sleeve height: 7mm Height incl. screw: 14mm	
	P35513	ZERAMEX®P6 Transfer Closed Tray, long SN	Sleeve height: 11mm Height incl. screw: 19mm	PEEK-CW30, Aluminum
	P35514	ZERAMEX®P6 Scanbody SN	Height: 10mm	




PRODUCT RANGE





ZERAMEX® Prosthetics - SN Platform

Auxiliaries	Art. No.	Name	Dimension	Material
	P35520	ZERAMEX®P6 Implant Replica SN	Height: 10mm	Aluminum
	P35550	ZERAMEX®P6 Provisional Screw SN	Length: 7.3mm	PEEK-CW30
	P35553	ZERAMEX®P6 Lab Screw SN	Length: 6.2mm	PEEK



ZERAMEX® Prosthetics - RN Platform

Soft Tissue Management	Art. No.	Name	Dimension	Material
	P36500	ZERAMEX®P6 Healing Cap RN	Height: 0.6mm	PEEK
	P36502	ZERAMEX®P6 Closure Screw RN	Height: 0.5mm	
	P36503	ZERAMEX®P6 Gingiva Former RN, 3 mm	Height: 3mm	
	P36504	ZERAMEX®P6 Gingiva Former RN, 4 mm	Height: 4mm	
	P36530	ZERAMEX®P6 Provisional RN (incl. screw)	AH: 7mm CH: 1mm	PEEK, PEEK-CW30

Impression Taking	Art. No.	Name	Dimension	Material
	P36510	ZERAMEX®P6 Transfer Open Tray RN	Sleeve height: 11mm Height incl. screw: 20mm	PEEK-CW30, Aluminum
	P36512	ZERAMEX®P6 Transfer Closed Tray RN	Sleeve height: 7mm Height incl. screw: 14mm	PEEK-CW30, Aluminum
	P36513	ZERAMEX®P6 Transfer Closed Tray, long RN	Sleeve height: 11mm Height incl. screw: 18mm	PEEK-CW30, Aluminum
	P36514	ZERAMEX®P6 Scanbody RN	Height: 10mm	PEEK, PEEK-CW30









Auxiliaries	Art. No.	Name	Dimension	Material
	P36520	ZERAMEX®P6 Implant Replica RN	Height: 10mm	Aluminum
	P36550	ZERAMEX®P6 Provisional Screw RN	Length: 6.1mm	PEEK-CW30
	P36553	ZERAMEX®P6 Lab Screw RN	Length: 5.2mm	PEEK
	P18550	ZERAMEX®P6 Abutment Planning Kit		PEEK

ZERAMEX® Tools














	Art. No.	Name	Material
	P35601	ZERAMEX®P Rose Drill ø 2.2 mm	Stainless steel
	P35602	ZERAMEX®P Pilot Drill ø 2.2 mm	
	P35633	ZERAMEX®P Drill ø 2.8 mm	
	P36633	ZERAMEX®P Drill ø 3.5 mm	
	P37633	ZERAMEX®P Drill ø 4.2 mm	
	P35622	ZERAMEX®P Profile Drill SN ø 2.8 mm	
	P36622	ZERAMEX®P Profile Drill ø 3.5 mm	
	P35650	ZERAMEX®P Depth Gauge ø 2.0 mm	
	P36650	ZERAMEX®P Depth Gauge ø 3.0 mm	
	P35620	ZERAMEX®P Screw Tap ø 3.3 mm	
	P36620	ZERAMEX®P Screw Tap ø 4.1 mm	
	P37620	ZERAMEX®P Screw Tap ø 4.8 mm	
	P38619	ZERAMEX® Prosthetic Wrench, 19 mm	
	P38623	ZERAMEX® Prosthetic Wrench, 23 mm	
	T38620	ZERAMEX® Pickup 20 mm	
	T38622	ZERAMEX® Rescue Pickup	
	T38625	ZERAMEX® Pickup 25 mm	
	T38626	ZERAMEX® Alignment Tool	
	C7650	Drill extension	




















PRODUCT RANGE

ZERAMEX® Tools

	Art. No.	Name	Material
	P35453	ZERAFIX™ Placeholder SN 3.0mm	PEEK
	P36453	ZERAFIX™ Placeholder RN 3.3mm	
	P38645S	ZERALOCK™ Abutment Holder Set SN & RN	PEEK
	P38644S	ZERAFIX™ Abutment Holder Set SN & RN	
	P48850	ZERAMEX®P Surgery Tray	Stainless steel
	P48854	ZERAMEX®P Surgery Tray, fully equipped with ratchet	
	P48870	ZERAMEX®P Surgery Tray "mini"	
	P48932	ZERAMEX®P Ratchet Adapter Unit Short	
	P48935	ZERAMEX®P Surgical Ratchet	
	P48860	ZERAMEX®P Prosthetics Tray	
	P48865	ZERAMEX®P Prosthetics Kit, incl. ratchet	

ZERAMEX® Locator® Accessories

	Art.No.	Name	Material
	T38219	LOCATOR® Processing Denture Cap (up to 10° divergence per implant, clear, pink, blue), 2pcs.	Ti/Nylon/LDPE
	T38220	LOCATOR® Processing Denture Cap (up to 20° divergence per implant, green, orange, red), 2pcs.	Ti/Nylon/LDPE
	T38222	LOCATOR® Blockout Spacer, 20pcs.	Silicone rubber
	T38224	LOCATOR® Impression Coping, 4pcs.	Al housing with LDPE lining
	T38226	LOCATOR® Replica Ø4mm, 4pcs.	Aluminum
	T38227	LOCATOR® Replica Ø5mm, 4pcs.	Aluminum
	T38228	LOCATOR® Core Tool	Stainless steel
	T38230	LOCATOR® Replacement clear, 4pcs. Retention force: 2.27kg	Nylon
	T38231	LOCATOR® Replacement pink, 4pcs. Retention force: 1.36kg	
	T38232	LOCATOR® Replacement blue, 4pcs. Retention force: 0.68kg	
	T38233	LOCATOR® Replacement green, 4pcs. Retention force: 1.81kg	
	T38234	LOCATOR® Replacement orange, 4pcs. Retention force: 0.91kg	
	T38235	LOCATOR® Replacement red, 4pcs. Retention force: 0.45kg	

	0050		Note (observe accompanying documents)
	Order number		Expiry date:
	Batch ID		Follow the instructions for use
	Do not reuse		Manufacturing date:
	Steam sterilized		Manufacturer
	Unsterile		Do not use if packaging is damaged
	Keep away from sunlight		The products of the portfolio are FDA approved
	Keep dry		Prescription use only
	Implant small ø3,3mm SN (Platform ø4,0mm)		Implant small ø4,1mm RN (Platform ø4,8mm)
	Implant small ø4,8mm RN (Platform ø4,8mm)		